Testsimplets®

Roche Diagnostics

(D)

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(Eng)

Intended use: Testsimplets® is an in vitro diagnostic rapid assay for staining and visualization of the differential blood count* under the microscope. Preparation of a blood film is no longer necessary¹².

Test principle: The test consists of a pre-stained slide. The stain contains two dyes, cresyl violet acetate and new methylene blue. The specific affinity of cell structures to the dyes results in different staining which permits classification of cells. The amounts of dye and the mixing ratio of the dyes in the pre-stained layer are constant and standardized, thus allowing reliable stains.

Contents of the pack: 50 ready-to-use slides pre-stained with 2.1 µg/cm² cresyl violet acetate and 1.0 µg/cm² new methylene blue.

Take the cassette out of the aluminium foil packing and open by lifting off the blue lid (Fig. 1). The slides can now be removed from the cassette by hand (Fig. 2) or with the aid of the dispenser (Fig. 3). **Be careful not to touch the pre-stained area of the slide.**

50 dust-free cover glasses (24 x 36 mm): To remove the cover glasses slit the label on the small cassette along the broken line and lift up the lid vertically in the direction of the arrow.

Storage and stability: + 2 °C to + 30 °C

Shelf-life: When kept in the original pack Testsimplets® are stable up to the date specified on the pack. Once the pack has been opened, high atmospheric humidity and large fluctuations in temperature may lead to sporadic formation of crystals in the pre-stained layer. Although this does not usually affect the performance of the test, such storage conditions should nevertheless be avoided.

Precautions and warnings: Correct evaluation of the blood count requires proper instruction and practice.

The cell layers will be too thick if too much blood has been applied and/or if dust has settled on the slide or cover glass. Cell differentiation in such cases is no longer possible.

Anticoagulants: EDTA is the recommended anticoagulant for venous blood. Venous blood containing EDTA must be kept at room temperature for no more than 3 hours before use (do not refrigerate)³.

Test procedure:

1. Apply a small drop of blood (about 3 µL) to the centre of a cover glass. 3 µL blood corresponds to a drop of blood that is about 3-4 mm in diameter on the slide. If the drop of blood is too large, the resulting layer of blood may be too thick and differentiation may not be possible.
- Capillary blood: Carefully collect a small drop of capillary blood from the fingertip with the cover glass.

Venous blood: Apply directly to the cover glass, for example using a 3-µL piston pipette.

2. Place the cover glass on the pre-stained area of the slide so that the drop of blood is in the centre of the area. The blood can also be carefully applied directly to the centre of the pre-stained area of the slide; the cover glass is then immediately placed on top.

3. If the blood does not spread out evenly and quickly, a sufficiently thin layer can be obtained by exerting mild pressure on the cover glass. Press down the centre of the cover glass with the point of a pencil and draw outwards in different directions.

4. After a 15-minute waiting period the preparation can be differentiated under oil immersion at 800x to 1000x. The prepared slide is stable for at least 4 hours at room temperature. If microscopic differentiation has to be delayed for 4-24 hours the slide should be refrigerated immediately after preparation.

Evaluation: For evaluation, select a site on the blood film where the cells are spread thinly enough to allow easy identification. Count 100 leukocytes, scanning the thin layer of the film in a meandering pattern.

While scanning the slide it is advisable to keep adjusting the focus of the leukocytes by turning the micrometer screw back and forth. In doubtful cases - such as may occur with the identification of monocytes, band neutrophils and basophils - it is particularly important to scan the entire cell, bringing individual cell regions into sharp focus.

Limitations of the method:

- The stability of the prepared film is limited to 4 hours at room temperature or 8 hours in the refrigerator.

- The method is of limited suitability for the red cell count.

- Testsimplets® do not allow documentation of the blood count.

Test results:

Accuracy: In the normal white blood count all forms of leukocytes are clearly identified and there is good agreement between the results of Testsimplets® and Pappenheim staining¹².

Reproducibility: For a normal white blood count within-series precision, day-to-day precision and laboratory-to-laboratory precision are the same for Testsimplets® and Pappenheim staining.

Interference: Use of citrated blood may lead to overstaining. Oxalate and fluoride are not suitable, likewise heparinate and blood samples from patients on heparin therapy.

* Further applications, including spermatozoal staining, urine cytology, cancer cytology and CSF cytology, are described in the scientific brochure published by Roche Diagnostics.

Presentation: 1 cassette containing 50 pre-stained slides and 1 cassette containing 50 dust-free cover glasses.

REF 191574

For an explanation of the symbols used and a list of references please refer to the end of this insert.

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Testsimplets is a trademark of Waldeck GmbH & Co. KG
Waldeck GmbH & Co. KG