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Manufacturer, classification and Intended Use**Producer and production location:**

Producer: Nipro Corporation, Japan.

Production place: Thailand (Nipro Thailand Corporation Limited).

European representative: Nipro Europe N.V. Belgium.**Classification:** Medical device Class IIa CE mark 0123 in accordance with the European Directive 93/42/EEC.**Intended use:** venipuncture and subsequent suction or injection of fluids from / into the vascular system.**Compliance:****Product:**

- EU Directive 93/42 / EEC;
- sterile hypodermic needles in accordance with ISO 7864;
- stainless steel needle cannula for the production of medical devices: ISO 9626;
- Designed in accordance with the harmonized standards to ensure a complete connection to other devices, provided that the other devices that they also comply with harmonized standards as EN 1707, ANSI / HIMA MD 70.1;
- The color coding for identification complies with ISO 6009;
- Medical device equipped with a safety mechanism against accidental punctures in compliance with the ISO23908 requirements;
- EN980 / ISO15223-1;
- complies Sterilization: ISO11135-1 / ISO11135 / EN556-1 / F.U. ;
- Non-toxic / Biocompatibility tested in accordance with ISO 10993;
- F.U. current ed.

Production:

- ISO 13485;
- Environmental Certification ISO14001;
- Environmental Policy compliant ISO14001-140021 / 4/5;
- Code of Conduct (Code of Ethics) according OHSAS 18001 and Legislative Decree no. 231/2001.

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» *The new model butterfly needle with protection **SafeTouch Nipro** ensures excellent venipuncture with a minimum resistance to penetration and maximum operator safety. The design combines safety with a simple practical use, intuitive and safe.*

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» *The stainless steel **needle**, is in ultra-thin wall, honed ultrasonic, internal and external siliconizing, special sharpening ultrasonic" Lancet Point "triple facet to ensure an atraumatic puncture. The needle is of such strength that even if bent 90° and returned to the original position does not lose their functional characteristics;*

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» **Protector Safe Touch:** after use, the needle body runs coaxially to the needle, covering it completely and irreversibly;

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» *The tube in flexible PVC (DEHP-free) medical grade, without memory, shows a higher biocompatibility and eliminates the risk of bottlenecks or breakages.*

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» *The spacious and flexible **wings** provide an excellent grip and are color coded for quick and unique identification of the needle size. The needle size and brand are also stamped on the same wings;*

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» **Sterilization and Validity:** sterilized with Ethylene Oxide, validity 60 months.

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Features of the mechanism SafeTouch Security:

- *Stop Security System: once activated, the protector device is irreversible and covers completely the needle, ensuring protection against accidental punctures to the operators and employees for product disposal. The needle is encapsulated inside of the protector device thanks to the special design of the product. Even when driven by shocks or forced, the needle does not come from the protector device. Even in case of system failure, deliberate or accidental, the needle remains encapsulated in the protector device;*

- *Non-removable security system: unlike many safety butterfly needles available on market, the security system of the needles Nipro SafeTouch is an integral part of the device, not a pre-assembled accessory. For this reason it can not be removed either before or after activation. - Can be activated with only one hand: the second hand is free to exert digital pressure on the area of the venipuncture after extraction of the needle, as per protocol in use;*

- *The activation of the safety mechanism is confirmed by a "click";*

- *The security device does not interfere with the use of the device: In order to provide maximum comfort to the patient before activating the safety device does not in any way limit the field of vision and action from the operator;*

- *The safety device is easy and intuitive, handy and practical;*

- *The design of the protector device allows use much akin to traditional nonprotected device.*

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when the needle is fully retracted is activated the locking system, with perceived "click" acoustically and mechanically.

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For more specifications and instructions for use see the instructions provided with the device.

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» **Packaging:** individual blister, collected in boxes of 50 pieces and packages of 500 pieces. On the single blister, carton of 50 pieces and the packaging contains all the information needed for rapid and unambiguous identification of the product, even in Italian or in the form of pictograms in compliance with EN 980, including: content description, quantity, needle size, lot number, expiration date, indication of sterile product and method sterilization, indication of disposable product, manufacturer's representative and the European CE mark 0123. To ensure a good stability over time, it is recommended to store products in dry and away from heat sources at a temperature comprised between 5 and 45 ° C.

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