

**SCHEDA TECNICA PRODOTTO / TECHNICAL DATA SHEET****PORTAVETRINI IN CARTONE
CARDBOARD SLIDE HOLDER****BSI204****BSI205**

Codice Code	Descrizione Description	Coperchio Lid	Divisori Dividers	Unità di vendita Sales unit
BSI099	1 posto / 1 place	x	x	50
BSI100	2 posti / 2 places	x	x	50
BSI200	32 posti / 32 places	x	-	5
BSI201	4 posti / 4 places	x	x	50
BSI202	9 posti / 9 places	-	-	10
BSI203	10 posti / 10 places	x	x	5
BSI204	20 posti / 20 places	x	x	5
BSI205	20 posti / 20 places	-	x	10
BSI206	24 posti / 24 places	x	-	5

Usò previsto

Adatto per la conservazione temporanea ed il trasferimento di vetrini.

Per usi generici.

Caratteristiche

Realizzato in cartone pesante, economico e riutilizzabile. Gli scomparti sono incassati per garantire che il campione non possa essere toccato quando il coperchio è chiuso.

Adatto per vetrini standard. (Dim. 26x76 mm; spessore 1.0±0.2 mm).

Precauzioni

Evitare ambienti umidi.

Non autoclavare.

Conservazione e durata di conservazione

Conservare in un ambiente pulito e asciutto a temperatura ambiente 15-30°C.

Intended use

Suitable for temporary storage and slide transfer.

For general uses.

Characteristics

Made of heavy cardboard, economical and reusable.

The compartments are recessed to ensure that the sample cannot be touched when the lid is closed.

Fits standard slides. (Dim. 26x76 mm; thickness 1.0±0.2 mm).

Precautions

Avoid humid environments.

Do not autoclave.

Storage and shelf life

Store in a clean, dry environment at room temperature 15-30°C.